

Indication Guide

Suitable for All Indications

Indications	Veneers	Inlays/onlays	Full contour anterior crown	Anterior bridge (3 units)	Full contour posterior crown	Posterior bridge (3 units)	Full contour bridge (4 units)	Full contour bridge (≥4 units)	Anterior coping	Anterior coping bridge (3 units)	Posterior coping	Reduced posterior bridge (3 units)	Reduced posterior bridge (4units)	Reduced posterior bridge (≥4units)
HT White	●	●	●	●	●	●	●	●	★	★	★	★	★	★
ST White	●	●	●	●	✓	✓	✓	★	★	★	★	★	★	★
ST Pre-shaded	●	●	●	●	✓	✓	✓	★	★	★	★	★	★	★
ST Multilayer	●	●	●	●	✓	✓	✓	★	★	★	★	★	★	★
TT ONE White	●	●	✓	✓	★	★	★	✓	✓	✓	✓	✓	✓	●
TT ONE Pre-shaded	●	●	✓	✓	★	★	★	✓	✓	✓	✓	✓	✓	●
TT ONE Multilayer	●	●	✓	✓	★	★	★	✓	✓	✓	✓	✓	✓	●
TT White	●	●	★	★	★	★	×	×	✓	✓	✓	✓	×	×
TT Multilayer	●	●	★	★	★	★	×	×	✓	✓	✓	✓	×	×
Explore Functional	●	●	✓	✓	★	★	★	✓	✓	✓	✓	✓	✓	✓
Explore Esthetic	●	●	★	★	★	★	●	×	✓	✓	✓	✓	●	×

★ optimum ✓ available ● available but not recommended × unavailable