

UPCERA HiFi Fire

Firing Furnace

Intuitive Operating System,
Customized Exclusive
Security Program



High-Quality
Craftsmanship



Durable
Stability

Applications	with vacuum	without vacuum
Fire	✓	✓
Press	-	-
Sinter	-	-
Glaze	✓	✓
Crystallize	✓	✓
Infiltrate	✓	✓

SIZE (W x H x D)	383 x 653 x 343 mm
WEIGHT	17.5 kg
POWER CONSUMPTION	95-105V / 50-60Hz or 110-125V / 50-60 Hz or 220-240V / 50-60 Hz
P (max)	1500 W
T (max)	1200 °C

- Even heat distribution through closed chamber design
- Automatic control of the drying and cooling to optimize aesthetics
- UPCERA HiFi Fire's specific firing chamber concept for excellent firing results in all temperature ranges
 - Thermocouple made of platinum/platinum-rhodium, purest quality with an accuracy of measurement far below 0.5°C
 - Individually selectable drying positions of the lift
 - Variable vacuum switch-on point, independent of preheating temperature
- Access to more than 100.000 firings, pressing, and sintering programs of all manufacturers of dental ceramics online

