


UPCERA

HiFi

FIRE

PRESS

SINTER



UPCERA

HIF

SIZE (BxHxT)	383 x 780 x 343 mm
WEIGHT	24 kg
POWER SUPPLY	220 – 240 V/50 – 60 Hz
P_{max}	2.000 W
T_{max}	1.560 °C

- ☪ Already prepared for ultrafast cycles below 1 hour
- ☪ DEKEMA patented high-temperature precision thermocouple design, suitable for speeding without overshooting temperature
- ☪ DEKEMA patented tray system with up to 2 levels (Ø 100 mm per level), invented for speeding or easy handling of up to 80 units
- ☪ No contamination and no discoloration
- ☪ Absolutely clean – no purge cycles necessary
- ☪ Nonstop – start subsequent cycles immediately
- ☪ DEKEMA AutoDry™
- ☪ Internet connectivity: profit from a secure backup of the process data, a quality management and fast service to minimize downtime