

UPCERA

HiFi

FIRE

PRESS

SINTER



UPCERA

HIF

SIZE (BxHxT)	383 x 653 x 343 mm
WEIGHT	22 kg
POWER SUPPLY	95 – 105 V/50–60 Hz or 110 – 125 V/50–60 Hz or 220 – 240 V/50–60 Hz
P_{max}	1.500 W
T_{max}	1.200 °C

- ☪ Automatic control of the drying and cooling to optimize aesthetics
- ☪ Multiple Press Time Automatic for optimal pressing results
- ☪ Besides all common ring systems, press up to 20 g with the exclusive DEKEMA trixpress™ ring system
- ☪ Press up to 3 different shades in one go, optimizing throughput
- ☪ Even heat distribution over all surfaces of the investment ring which provides for exceptional aesthetics
- ☪ DEKEMA AutoDry™
- ☪ Internet connectivity: profit from a secure backup of the process data, a quality management and fast service to minimize downtime